ES

IT

DE

EN

#### Descripción:

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#### Description:

**Alpha-dent®** Light Cure Hybrid Composite Restorative Material is a Bis-GMA polymeric matrix and a 56% by volume blend of inorganic radiopaque fillers with particle range from 0.02 to 4 microns. **Alpha-dent®** Light Cure Composite is available in 9 Vita\* shades plus one opaquer and an incisal shade.

#### Indications for Use:

**Alpha-dent®** Light Cure Composite Restorative Material is recommended for Class III and Class V restorations and for limited use in Class I restorations in pre-molars and selected Class IV restorations where esthetics is of primary importance.

#### Contraindication:

**Alpha-dent®** Light Cure Composite Restorative Material is contraindicated in patients with a history of severe allergic reaction to methacrylate resins.

#### Warnings:

- Avoid contact of **Alpha-dent®** etchant and **Alpha-Etch37®** with soft tissue, eyes and skin. Both contain phosphoric acid, which is corrosive. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek immediate medical attention.
- Do not use in the presence of eugenol based materials as complete setting may not occur.

#### General Precautions:

- U. S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.

#### Adverse Reaction:


Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.

#### Instructions for Use:

- Prepare the tooth. Isolate and dry the tooth thoroughly.
- If using **Alpha-dent®** etchant or **Alpha-Etch37®** etch enamel for 15 seconds or if applied to dentin etch for 10 seconds; never leave etchant on tooth for more than 60 seconds. Lightly agitating the etchant during its application can enhance its action.
- Rinse thoroughly with water and dry tooth completely.
- If using **Alpha-dent®** Light Cure Bonding Adhesive, thoroughly dry the etched tooth. Dispense two drops of **Alpha-dent®** Light Cure Bonding Adhesive onto a dispensing pad or other suitable surface. Thoroughly saturate a brush tip and generously apply a uniform layer onto the prepared tooth. Light cure the bonding adhesive by applying a dental halogen light for 15-20 seconds. The surface should appear glossy. If not, repeat application.
- Select the appropriate shade of **Alpha-dent®** Light Cure Composite. Apply the composite directly onto the prepared tooth in increments no larger than 2 mm thickness.
- Using a visible light-curing unit, with a minimum output of 600mW/cm<sup>2</sup>, cure each incremental layer for at least 20 seconds.
- Finish and polish using standard techniques. A composite surface sealant can also be used to improve surface smoothness and seal micro cracks and surface imperfections caused by the finishing procedures.

#### Storage:

- Alpha-dent®** Light Cure Composite is formulated to be used at room temperature. See outer package for expiry date. Do not freeze. Recommended storage conditions are between 2°C/36°F to 27°C/80°F, bringing to room temperature before using (23°C/73°F).
- Do expose restorative materials to elevated temperatures or intense light.
- If using unit-dose tips, which are intended for single patient use, discard after each patient. If using syringes, always replace cap after each use.

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