Descrpción:

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Description:

Alpha-Bond® Single Component Light Cure Bonding Adhesive with Dual Cure Activator (purchased separately) is formulated to bond to multiple substrates for both direct and indirect bonding procedures. It contains the hydrophilic monomer PMGDM that allows penetration into the tubules, intertubules and peritubular dentin. Our bonding adhesive offers you the flexibility and confidence you need in a bonding adhesive. For a single step light cure adhesive, use the bottle marked bonding adhesive. By combining the dual cure activator with the adhesive you can confidently bond to self-cure composites and areas inaccessible to light.

Indications for use:

- Direct Bonding of Composite Resins and Compomers
- Veneers
- Composite, Ceramic and Metal repairs
- Indirect restorations: inlays, onlays, crown and bridge retainers
- Endodontic post cementation
- Direct Bonding of Composites and Core Build-up Materials

Contraindications:

Alpha-Bond® is contraindicated in patients with a history of allergic reaction to methacrylate resins.

Warning:

- **Alpha-Bond®** is highly flammable, keep away from sources of ignition.

General Precautions:

- U. S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed to light.

- Always replace cap after each use.

Adverse Reaction:

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.

Instructions for use Light Cured Materials

1. Prepare the tooth. Isolate and dry tooth with oil-free compressed air.
2. If using **Alpha-dent®** etchant or **Alpha-Etch37®**, apply etchant to the prepared enamel for 15 seconds, or if applied to dentin etch for 10 seconds. Lightly agitating the etchant during its application can enhance its action. Never leave the etchant on the tooth for more than 60 seconds.
3. Rinse thoroughly with water. Remove excess water with suction followed by blotting with suitable dry lint-free material, leaving the dentin and enamel with a non-shiny, matte surface yet moist.
4. Dispense two drops of **Alpha-Bond®** onto a suitable surface. Thoroughly saturate an applicator tip and generously apply two coats onto the tooth, remove excess adhesive by suctioning, being careful to keep a continuous coating on the tooth.
5. Gently apply air for at least 5 seconds to evaporate the solvent, making sure you do not apply too much air to thin the adhesive. When properly dried the adhesive layer will not move with air and the surface should appear shiny. If the surface is not shiny it may not be properly sealed and additional applications are recommended.
6. Using a light-curing unit (minimum output of 600 mW/cm2) cure the material for at least 10 seconds. The light-curing unit tip should be placed in close proximity to surface being treated.
7. Immediately place restorative material following the manufacturer's instructions.

Veneers:

1. Follow steps 1-6 from Direct Bonding Composite Resins and Compomers to Dentin/Enamel.
2. Treat bonding surface of porcelain or ceramic veneer according to manufacturer's or laboratory's instructions (i.e. etching, mechanical roughening, silanating, etc.)
3. Apply luting cement to the veneer, seat the veneer, and cure the luting cement according to the manufacturer's instructions.

Composite, Ceramic and Metal Repairs:

1. Prepare and roughen the fractured surface with a diamond bur or by sandblasting. Create mechanical retention where possible. Rinse, isolate, and dry thoroughly with oil-free compressed air. If repairing composite and/or metal proceed to step 3.
2. If repair is to porcelain, remove all loose pieces and contour porcelain edge. Isolate area with rubber dam and etch surface with hydrofluoric acid according to manufacturer's instructions. Rinse thoroughly with water and dry with air. Apply silane to etched porcelain surface.
3. Follow steps 1-6 from Direct Bonding Composite Resins and Compomers to Dentin/Enamel.
4. Place the restorative material over the cured **Alpha-Bond®** and cure the luting cement according to the manufacturer's instructions.

Instructions for use Self and Dual Cure Materials - Mixing Ratio - Two drops of Bonding Adhesive to One drop of Activator:

1. Indirect restorations: inlays, onlays, crown and bridge retainers:
 1. Treat bonding surface of porcelain or ceramic restoration according to manufacturer's or laboratory's instructions (i.e. etching, mechanical roughening, silanating, etc.)
 2. Follow steps 1-3 from Direct Bonding Composite Resins and Compomers to Dentin/Enamel.
 3. Dispense two drops of **Alpha-Bond®** Adhesive and one drop of **Alpha-Bond®** Activator onto a suitable surface. Mix contents for 1-2 seconds with a clean applicator tip. Thoroughly saturate an applicator tip and generously apply two coats onto the tooth and bonding surface of the porcelain/ceramic restoration, remove excess adhesive by suctioning being careful to keep a continuous coating on the surface. This should be followed by gently applying air for at least 5 seconds to evaporate the solvent, making sure you do not apply too much air to thin the adhesive. When properly dried the adhesive layer will not move with air and the surface should appear shiny. If the surface is not shiny it may not be properly sealed and additional applications are recommended. Do not light cure the **Alpha-Bond®** Dual Cure mixture.
 4. Prepare and apply cement according to the manufacturer's instructions.

Endodontic Post Cementation:

1. Treat bonding surface of post according to manufacturer's or laboratory's instructions (i.e. etching, mechanical roughening, silanating, etc).
2. Prepare the post space, isolate and dry thoroughly using air, suction or paper points.
3. Etch the post space using the **Alpha-Etch37®** technique. Apply etchant to the prepared post space. Etch dentin for 10 seconds; never leave etchant on tooth for more than 60 seconds.
4. Rinse thoroughly with water. Remove excess water in post space with a brief burst of air, suction or paper points, leaving the dentin with a non-shiny, matte surface yet moist.
5. Dispense two drops of **Alpha-Bond®** Adhesive and one drop of **Alpha-Bond®** Activator onto a suitable surface. Mix contents for 1-2 seconds with a clean applicator tip. Thoroughly saturate an applicator tip or paper point and generously apply two coats into the preparation, remove excess adhesive by suctioning being careful to keep a continuous coating on the tooth. This should be followed by gently applying air for at least 5 seconds to evaporate the solvent, making sure you do not apply too much air to thin the adhesive. When properly dried the adhesive layer will not move with air and the surface should appear shiny. If the surface is not shiny it may not be properly sealed and additional applications are recommended. Do not light cure the **Alpha-Bond®** Dual Cure mixture.
6. Apply a single coat of **Alpha-Bond®** Dual Cure Mixture to the post with the same applicator tip. Evaporate the solvent with a gentle stream of air for 5 seconds, making sure you do not apply too much air to thin the adhesive. If the surface is not shiny on the post reapply a second coat of **Alpha-Bond®** Dual Cure Mixture and evaporate the solvent with a gentle steam of air for 5 seconds. Do not light cure the **Alpha-Bond®** Dual Cure mixture.
7. Prepare and apply cement to surface of post and/or into the post preparation according to the manufacturer's instructions.

Direct Bonding of Self Cure/Dual Cure Composites and Core Build-up Material:

1. Follow steps 1-3 from Direct Bonding Composite Resins and Compomers to Dentin/Enamel.
2. Dispense two drops of **Alpha-Bond®** Adhesive and one drop of **Alpha-Bond®** Activator onto a suitable surface. Mix contents for 1-2 seconds with a clean applicator tip. Thoroughly saturate an applicator tip and generously apply two coats of

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