



alpha-dent®

# ORANGE SOLVENT

WITH LANOLIN

8 fl. oz.

REF 005.0102.001

**USES:** For the removal of impression materials from hands, face, mixing slabs and instruments.

**DIRECTIONS:** Moisten cloth or swab with ORANGE SOLVENT and wipe away impression material.

**INGREDIENTS:** Mineral Oil, Limonene, Lanolin.



**DENTAL  
TECHNOLOGIES**

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA  
[www.dentaltech.com](http://www.dentaltech.com)

PM-2419/REV001