



TECHNICAL DATA SHEET

Product: Alpha-Dent® Alginate Flavorings

Product Description: Highly concentrated liquid flavorings that add a unique and pleasant flavor to alginate impression materials.

Indications for Use: Adding flavor to alginate impression materials.

Appearance: Clear colorless to dark-colored liquid.

Curing Mechanism: Not applicable.

Shelf-Life: 30 months

Filler Type: None

Main Composition (%w/w):

| | |
|--------------------|---------|
| Solvent | 80-90 |
| Sodium Saccharin | <5 |
| Surfactant | <1 |
| Flavorings: | 10 – 20 |
| - Propylene Glycol | |
| - Alcohol | |
| - Water | |
| - N&A Flavors | |

Physical /Mechanical Properties:

| | |
|-------------------|------------------------------------|
| Appearance: | Medium viscosity liquid |
| Color: | Characteristic of the given flavor |
| Flavor: | Characteristic of the given flavor |
| Specific gravity: | 0.99 to 1.09 |

Regulatory Classification:

This product is in conformity with the following standard(s) or other normative documents:
European Standard BS EN 1641:2009 – Medical Devices for Dentistry – Materials
21 CFR 820 – Medical Device Quality System Regulation

Signature/Date: John Pontikis / 10/18/18

Name: John Pontikis
Position: Management Representative