



TECHNICAL DATA SHEET

Product: Alpha-Ease® Desensitizer

Product Description: A one-step product for treating and preventing dentinal hypersensitivity. It works quickly and effectively to seal the dentinal tubules.

Indications for Use: Under any indirect resin-bonded inlay, onlay, crown or veneer. Under amalgam and composite restorations. In conjunction with other dentin bonding adhesives. For treating hypersensitive cervical areas and reducing/eliminating postoperative sensitivity.

Appearance: Clear colorless to pale yellow liquid

Shelf-Life: 30 months

Main Composition (%w/w):

Water	40-49
Hydroxyethylmethacrylate	20-50
Glutaraldehyde	<6
Sodium Fluoride	<1

Physical /Mechanical Properties: (typical values)

% fluoride ion	0.18
Refractive Index	1.39
Specific Gravity	1.04

Regulatory Classification:

This product is in conformity with the following standard(s) or other normative documents:

European Standard EN 1641:2009 – Medical Devices for Dentistry – Materials
Code of Federal Regulations CFR 21 Part 820 – Medical Device Quality System Regulation

Signature/ Date:  / 10/23/18

Name: John Pontikis
Position: Management Representative