



# TECHNICAL DATA SHEET

**Product:** Alpha-Dent® Orange Solvent

**Product Description:** Alpha-Dent® Orange Solvent is a mild orange-based solvent used for the cleaning and removal of cements and pastes from the patient's skin and tissue, as well as cleaning mixing slabs and instruments.

**Indications for Use:** Designed to be used for cleaning and removal of cements and pastes from the patient's skin and tissue, as well as for cleaning mixing slabs and instruments.

**Appearance:** Clear, light yellow liquid.

**Curing Mechanism:** Not applicable.

**Shelf-Life:** 48 months

**Filler Type:** None

**Main Composition (%w/w):**

Natural and artificial oils	100
-----------------------------	-----

**Physical Properties: (Typical Values)**

Appearance	clear liquid
Color	light yellow
Odor	lemon
Specific Gravity	0.84
Refractive index	1.467

**Regulatory Classification:**

This product is in conformity with the following standards or other normative documents:

Code of Federal Regulations 21 CFR §701.3(a) – Cosmetic Labeling

**Signature/ Date:** John Pontikis / 10/25/18

**Name:** John Pontikis  
**Position:** Management Representative