



TECHNICAL DATA SHEET

Product: Alpha-Pro® White Varnish

Product Description: Alpha-Pro® White Varnish contains 5% sodium fluoride. It is used in the treatment of dentinal and post operative sensitivity. Alpha-Pro® White Varnish, when applied, practically becomes invisible, unlike traditional varnishes that cause unsightly yellow coloring of the teeth. Alpha-Pro® White Varnish sets immediately on contact with saliva, leaving a smooth, thin layer. It is available in bubblegum, melon, mint, and strawberry flavors.

Indications for Use: For treatment of dentinal and post operative sensitivity.

Appearance: White to off-white, medium viscosity liquid.

Curing Mechanism: N/A

Shelf-Life: 30 months

Filler Type: N/A

Main Composition (%w/w):

Colophony:	45 - 60
Ethyl Alcohol:	20 - 35
Sodium Fluoride:	5
Emulsifier:	2 - 7
Flavor:	2 - 5
Xylitol:	<1
Sweetener:	<1

Physical /Mechanical Properties:

Appearance:	Medium viscosity liquid
Color:	White to off-white
Odor:	Characteristic of ethanol and specific flavor
Flavor:	Bubblegum, melon, mint or strawberry
% Sodium Fluoride:	5

Regulatory Classification:

This product is in conformity with the following standard(s) or other normative documents:

Code of Federal Regulations CFR 21 Part 820 – Medical Device Quality System Regulation;
Code of Federal Regulations CFR 21 Part 872.3260 – Cavity varnish.

Signature/ Date: John Pontikis / 10/25/18

Name: John Pontikis
Position: Management Representative