



# ALPHA-DENT® LIGHT CURE BONDING ADHESIVE

ALPHA-DENT® PROVIDES ↓

- BONDS TO DENTIN & ENAMEL
- IDEAL VISCOSITY
- COST EFFECTIVENESS
- AVAILABLE WITH ETCH LIQUID OR GEL
- GOOD BOND STRENGTH

## LIGHT CURE BONDING ADHESIVE



### LIGHT CURE BONDING ADHESIVE

Alpha-Dent® Light Cure Bonding Adhesive is designed for bonding composite resin materials to dentin and enamel. The ideal viscosity allows the material to penetrate into microtubules to enhance mechanical adhesion.

### Physical/ Mechanical Properties

Depth of Cure, mm (actual)	6.0
Bond Strength to Dentin, MPa	15
Bond Strength to Enamel, MPa	19.9



#### ALPHA-DENT®

REF 206.0000.001 ETCHANT LIQUID KIT

- KIT CONTAINS:
- 2 x 3mL BONDING ADHESIVE
- 7mL ETCH LIQUID
- DISPENSING PAD
- BRUSHES

REF 206.0001.001 ETCHANT GEL KIT

- KIT CONTAINS:
- 2 x 3mL BONDING ADHESIVE
- 2 x 4g ETCH GEL SYRINGES
- DISPENSING PAD
- BRUSHES

### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX