



# ALPHA-FLOW® FLOWABLE COMPOSITE

ALPHA-FLOW® PROVIDES ↓

- CONTAINS FLUORESCENCE
- EXCELLENT MECHANICAL PROPERTIES
- RADIOPAQUE
- IDEAL VISCOSITY
- CONTAINS FLUORIDE



## FLOWABLE COMPOSITE

**Alpha-Flow®** is a radiopaque fluoride containing light cure resin-based microhybrid flowable composite. The specially formulated material contains 65% by weight (38% by volume) inorganic filler. It incorporates both a nano-sized filler for increased flexibility and superior polishability and 0.7 micron filler for excellent strength and wear resistance. This special combination provides a material with extreme versatility for use in a wide range of procedures, including Class III and V cavity figurations and small Class I and II restorations. **Alpha-Flow®** contains fluorescence to provide duplication of natural tooth surface.

Available in six popular Vita® shades (A1, A2, A3, A3.5, B1, B2) and Universal Opaquer

Compressive Strength	206 MPa
Flexural Strength	110 MPa
Diametral Tensile Strength	42 MPa
Filler Particle Range ( $\mu$ )	0.02 - 0.70
Average Filler Particle Size ( $\mu$ )	0.65

<p><b>ALPHA-FLOW®</b></p> <p>REF 103.0003.002 STANDARD KIT</p> <p>REF 103.0004.002 CUSTOMKIT</p> <p>KIT CONTAINS:</p> <ul style="list-style-type: none"> <li>■ 4 X 2g SYRINGE</li> <li>■ 20 DISPENSING TIPS</li> </ul>
--

## DENTAL TECHNOLOGIES

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US AND CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX

