



ALPHA-DENT® LIGHT CURE BONDING ADHESIVE

ALPHA-DENT® PROVIDES ↓

- BONDS TO DENTIN & ENAMEL
- IDEAL VISCOSITY
- COST EFFECTIVENESS
- AVAILABLE WITH ETCH LIQUID OR GEL
- GOOD BOND STRENGTH

LIGHT CURE BONDING ADHESIVE



LIGHT CURE BONDING ADHESIVE

Alpha-Dent® Light Cure Bonding Adhesive is designed for bonding composite resin materials to dentin and enamel. The ideal viscosity allows the material to penetrate into microtubules to enhance mechanical adhesion.

Physical/ Mechanical Properties

Depth of Cure, mm (actual)	6.0
Bond Strength to Dentin, MPa	15
Bond Strength to Enamel, MPa	19.9



ALPHA-DENT®
REF 206.0000.001 ETCHANT LIQUID KIT
<ul style="list-style-type: none"> ■ KIT CONTAINS: ■ 2 x 3mL BONDING ADHESIVE ■ 7mL ETCH LIQUID ■ DISPENSING PAD ■ BRUSHES
REF 206.0001.001 ETCHANT GEL KIT
<ul style="list-style-type: none"> ■ KIT CONTAINS: ■ 2 x 3mL BONDING ADHESIVE ■ 2 x 4g ETCH GEL SYRINGES ■ DISPENSING PAD ■ BRUSHES

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX