



# ALPHA-GLAZE®

## LIGHT CURE SURFACE SEALANT

ALPHA-GLAZE® PROVIDES ↓

- WEAR RESISTANT FINISH
- CLEAR IN COLOR
- COMPATIBLE WITH MOST RESIN BASED COMPOSITES
- IDEAL VISCOSITY



### LIGHT CURE SURFACE SEALANT

**Alpha Glaze®** Surface Sealant contains light cure unfilled resin with low viscosity to be used as a finishing sealant for composites and resin cements. It is formulated to penetrate and seal the micro cracks and surface imperfections that are created by the finishing procedures. It can also be used to treat the enamel surface prior to applying Orthodontic Appliances.

### Physical/ Mechanical Properties

Specific Gravity	1.10-1.13
Depth to Cure, mm	6.0 (min)

#### ALPHA-GLAZE®

REF CODE	PRODUCT DESCRIPTION
403.0101.002	7ml Surface Sealant

#### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
 1.847.677.5500 PHONE  
 1.847.677.5502 FAX

