



ALPHA-DENT® SELF CURE BONDING ADHESIVE

SELF CURE BONDING ADHESIVE



ALPHA-DENT® SELF CURE BONDING ADHESIVE PROVIDES

BONDS TO DENTIN & ENAMEL

IDEAL VISCOSITY

COST EFFECTIVENESS

AVAILABLE WITH ETCH LIQUID OR GEL

4 YEAR SHELF LIFE

SELF CURE BONDING ADHESIVE

Alpha-Dent® Self Cure Bonding Adhesive is designed for bonding composite resin materials to dentin and enamel. The ideal viscosity allows the material to penetrate into microtubules to enhance mechanical adhesion.

Physical/ Mechanical Properties

Working Time at 23 °C, seconds	35
Setting Time at 37 °C, seconds	60
Bond Strength to Dentin, MPa	15
Bond Strength to Enamel, MPa	25



ALPHA-DENT®

REF 200.0000.001 ETCHANT LIQUID KIT

KIT CONTAINS:

- 3mL BASE ADHESIVE
- 3mL CATALYST ADHESIVE
- 7mL ETCH LIQUID
- MIXING PAD
- BRUSHES

REF 200.0001.001 ETCHANT GEL KIT

KIT CONTAINS:

- 3mL BASE ADHESIVE
- 3mL CATALYST ADHESIVE
- 2 x 4g ETCH GEL SYRINGES
- MIXING PAD
- BRUSHES

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]

1.847.677.5500 PHONE

1.847.677.5502 FAX