



ALPHA-DENT® CONDITIONERS

Metal, Porcelain, and Plastic Conditioners

ALPHA-DENT® CONDITIONERS PROVIDES



MULTIPLE SURFACES

IMPROVED BONDING

PREPARES SURFACES

CONVENIENT BRUSH ON CAP



METAL CONDITIONER

Alpha Dent® Metal Conditioner is to be used as a primer for amalgam, gold, or stainless steel. The conditioned metal surface will provide the optimal bonding surface for the adhesive.

PORCELAIN CONDITIONER

Alpha Dent® Porcelain Conditioner permits bonding versatility with patients that have bridgework or crowns. This conditioner is a silane coupling agent, which helps to chemically bond the adhesive to porcelain.

PLASTIC CONDITIONER

Alpha Dent® Plastic Conditioner prepares plastic brackets and buttons for direct bonding. Application of conditioner prior to adhesive application will provide an enhanced bond between the bracket base and the adhesive.



ALPHA DENT® CONDITIONERS

REF Code	Product
534.0101.001	15mL Metal Conditioner
542.0101.001	15mL Porcelain Conditioner
544.0101.001	15mL Plastic Conditioner

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]

1.847.677.5500 PHONE

1.847.677.5502 FAX