



ALPHA-DENT® POLYCARBOXYLATE CEMENT LUTING CEMENT

ALPHA-DENT® POLYCARBOXYLATE CEMENT PROVIDES ↓

- LUTING CEMENT
- SUPERIOR STRENGTH
- RESISTANCE TO DISSOLUTION
- ORAL RETENTION



POLYCARBOXYLATE LUTING CEMENT

Alpha Dent® Polycarboxylate Cement is specifically formulated for superior strength, resistance to dissolution, and oral retention. Our Luting Cement is impervious to oral fluids, and thus permanent under oral conditions. The material is recommended for cementing crowns, inlays, orthodontic bands, temporary fillings, fixed bridges, or as a base under any filling material.

Physical/ Mechanical Properties

Film Thickness	<25 µm
Net Setting Time, at 37 °C	2.5 - 3.0
Compressive Strength	70 Mpa, minimum
Acid-soluble Lead	< 100 mg/kg



POLYCARBOXYLATE CEMENT	
REF CODE	PRODUCT DESCRIPTION
306.0000.001	32g Powder, 16.5mL Liquid, Measuring Scoop, Mixing Pad

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX