



ALPHA-DENT® LIGHT CURE PIT & FISSURE

LIGHT CURE RESIN BASED PIT & FISSURE SEALANT

ALPHA-DENT® PROVIDES ↓

- WEAR RESISTANCE ■
- UNFILLED, RESIN BASED ■
- IDEAL VISCOSITY ■
- COST EFFECTIVENESS ■
- OPAQUE SHADE ■



LIGHT CURE RESIN BASED PIT & FISSURE SEALANT

The hard-setting wear resistant unfilled, resin based material has an ideal flowable viscosity for maximum coverage and protection to seal pits and fissures of teeth susceptible to caries. It's in an opaque shade to allow for easy identification during patient examinations.

Physical/ Mechanical Properties

Appearance	White Opaque Liquid
Depth of Cure mm/40 sec.	1.6



ALPHA-DENT® LC PIT & FISSURE

REF401.0101.001 LIGHT CURE SEALANT

- KIT CONTAINS
- 3mL LIGHT CURE OPAQUE SEALANT (WHITE)

REF401.0000.001 ETCHANT LIQUID KIT

- KIT CONTAINS
- 2 x 3mL SEALANT
 - 7mL ETCH LIQUID
 - DISPENSING PAD

REF401.0003.001 ETCHANT GEL KIT

- KIT CONTAINS:
- 2 x 3mL SEALANT
 - 2 x 4g ETCH GEL SYRINGES
 - DISPENSING PAD

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX