



ALPHA-CAL[®]

CALCIUM HYDROXIDE PASTE

ALPHA-CAL[®] PROVIDES ↓

- RADIOPAQUE ■
- pH > 12.0 ■
- SYRINGE READY ■
- DISPOSABLE TIPS ■



ALPHA-CAL CALCIUM HYDROXIDE PASTE

Alpha-Cal[®] is a calcium hydroxide paste with Barium Sulfate used as a root canal sealer for the temporary filling of root canals after conventional endodontic treatments. Alpha-Cal[®] is an aqueous radiopaque, high pH (> 12.0) paste packaged in a syringe ready to be dispensed through disposable tips into the root canal.



Physical/ Mechanical Properties

pH (10% aqueous dispersion)	12.4
Tensile Strength (thru # 20-gauge needle), inch	5
Specific Gravity	1.4

ALPHA-CAL [®]	
REF CODE	PRODUCT DESCRIPTION
414.0000.002	4 x 1g Syringe, 8 Syringe Tips

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX