



ALPHA-DENT® TRAY ADHESIVE

ALPHA-DENT® TRAY ADHESIVE PROVIDES



ACCURATE IMPRESSIONS

CONVENIENT BRUSH IN CAP

EASY TO USE



ALPHA-DENT® TRAY ADHESIVE

Alpha Dent® Tray Adhesive is used to obtain the most accurate impressions possible by adhering impression material to the tray. Since impression materials do not have adhesive properties, it is advisable to use **Alpha Dent® Tray Adhesive** on all trays before taking an impression.



Physical Properties: (Typical Values)

Appearance	viscous, slightly turbid liquid
Color	blue
Specific Gravity	0.86
Refractive index	1.438

ALPHA-DENT® TRAY ADHESIVE

REF Code	Product Description
540.0103.001	2 x 14ml (0.5 fl oz) Tray Adhesive

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX