




# ALPHA-DAM®

## LIGHT CURE GINGIVAL DAM

ALPHA-DAM® PROVIDES ↓

- CURES QUICKLY & EVENLY
- LOW EXOTHERMS
- APPLIES DIRECTLY
- REMOVES EASILY
- PROTECTS GINGIVAL TISSUE FROM CHEMICAL AGENTS



### LIGHT CURE GINGIVAL DAM

**Alpha-Dam®** serves as a barrier for in-office whitening treatments. The material cures quickly and evenly, while maintaining low exotherms (temperature rise of less than 12°C) to ensure patient comfort.

### Physical / Mechanical Properties (Typical Values):

Depth of Cure, mm/10sec. 7



ALPHA-DAM®	
REF CODE	DESCRIPTION
582.0101.002	3.5g Syringe with 5 Dispensing Tips

### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX