



# ALPHA-DENT® LIGHT CURE COMPOSITE

## LIGHT CURE HYBRID COMPOSITE

ALPHA-DENT® PROVIDES ↓

COST EFFECTIVENESS

CONTAINS FLUORESCENCE

RADIOPAQUE

NON-SLUMPING

FOURTEEN SHADES



### LIGHT-CURE HYBRID COMPOSITE

Alpha-dent® Light Cure Composite restorative is radiopaque. It has the versatility of the Hybrid composite while measuring high compressive strength with excellent polishability. This composite is appropriate for almost every type of restoration. The material is available with the Exact Match System, which uses a specially designed computer program to match the composite color and translucency to the VITA®<sup>1</sup> shade guide.

Compressive Strength	225 MPa
Flexural Strength	95 MPa
Diametral Tensile Strength	47 Mpa

#### ALPHA-DENT® LIGHT CURE COMPOSITE

REF CODE	PRODUCT DESCRIPTION
101.0001.001	(Custom) 7 x 4.5g Composite Syringes, 3ml Light Cure Bonding Resin, 7 ml Etchant Liquid, Brushes, Brush Holder, Dispensing Pad (Choice of shades)
101.0003.001	(Custom) 7 x 4.5g Composite Syringes, 7ml Light Cure Bonding Resin, 2 x 4g Etchant Gel Syringes & Tips, Brushes, Brush Holder, Dispensing Pad (Choice of shades)
101.0008.001	(Custom) 2 x 4.5g Composite Syringes, 3ml Light Cure Bonding Resin, 3ml Etchant Liquid, Brushes, Dispensing Pad (Choice of shades)
101.0010.001	(Custom) 4 x 4.5g Composite Syringes, 3ml Light Cure Bonding Resin, 7ml Etchant Liquid, Brushes, Dispensing Pad (Choice of shades)
101.0016.001	4 syringe with Etchant Gel

#### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX  
www.dentaltech.com

++ Choice of shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D4, Universal Opaquer, Incisal