




ALPHA-DENT® ALGINATE IMPRESSION MATERIAL

ALPHA DENT® ALGINATE IMPRESSION MATERIAL PROVIDES ↓

- HEAVY BODY
- EASY TO MIX
- SCOOP AND VIAL
- NORMAL OR FAST SET



ALGINATE IMPRESSION MATERIAL

Alpha-Dent® Alginate Impression Material is available in a fresh mint flavor, either in a normal or fast setting time. This heavy body material is designed to produce consistent impressions. Each can includes powder scoop and vial.

Physical/ Mechanical Properties

	FAST SET	NORMAL SET
Initial setting time (sec.)	90	155
Total working time (mm)	0.04	0.04
Recovery from deformation %	98	98
Strain in compression %	14.5	14.0
Compressive Strength Mpa	0.80	0.85
Compatibility to gypsum	Passes test	Passes test
Reproduction of details	Passes test	Passes test
Deterioration Mpa	0.6	0.5

ALPHA-DENT® ALGINATE IMPRESSION MATERIAL

REF CODE	PRODUCT DESCRIPTION
603.0103.001	Normal Set Mint 454g(1lb)
604.0103.001	Fast Set Mint 454(1lb)

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX