



ALPHA-DENT® ALGINATE FLAVORINGS


ALPHA DENT® ALGINATE FLAVORINGS PROVIDES

FUN FLAVORS

IMPROVES PATIENT'S EXPERIENCE

SIMPLE TO USE

GOOD FOR OVER 100 IMPRESSIONS



ALGINATE FLAVORINGS

Alpha-Dent® Alginate Flavorings adds a unique and pleasant flavor to each alginate mixture. The 19 fun flavors provide patients with a wide variety of choices, making the impression taking experience more enjoyable. Simple to use: just add 5 to 7 drops for every scoop of powder. The bottle contains 2 fl. oz. (59.2mL)

Physical/ Mechanical Properties

Appearance
Color
Flavor
Specific Gravity

Medium viscosity liquid
Characteristics of the given flavor
Characteristics of the given flavor
0.99 to 1.09

ALPHA-DENT® ALGINATE FLAVORINGS

REF CODE	PRODUCT DESCRIPTION
620.0102.001	Banana
620.0202.001	Blueberry
620.0302.001	Bubble Gum
620.0502.001	Cherry
620.0602.001	Cola
620.0902.001	Fruit Punch
620.1002.001	Grape
620.1202.001	Orange
620.1302.001	Peach
620.1402.001	Peppermint
620.1502.001	Pina Colada
620.1702.001	Raspberry
620.1802.001	Rootbeer
620.1902.001	Strawberry
620.2002.001	Vanilla
620.2102.001	Watermelon
620.2202.001	Apple
620.2802.001	Cotton Candy
620.3102.001	Cookie Dough

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX